



MICELLULA DNA EMULSION & PURIFICATION KIT



QUICK MANUAL

KIT VERSION 1.0, JULY 2010. SEE ALSO DETAILED VERSION OF THIS PROTOCOL.

A

EMULSION REACTION SET-UP

1

1. Emulsion - Setting up the Oil Phase

- Prepare 300 μ l Oil Surfactant Mixture per reaction (50 μ l water phase). Mix
220 μ l Emulsion Component 1 (73 % [v/v])
20 μ l Emulsion Component 2 (7 % [v/v])
60 μ l Emulsion Component 3 (20 % [v/v])
- Mix thoroughly by vortexing.
- Keep on wet ice until further usage.



2

2. Emulsion - Setting up the Water Phase

- Prepare 50 μ l Water Phase per Emulsion reaction (e.g. ePCR).
Mix template DNA and enzyme without BSA according to application requirements.
- Remove aliquot for "open" (unemulsified) reaction control (optional).
- Add BSA. Required amount: -----

3

3. Create Emulsion Reactions

- Mix 300 μ l Oil Surfactant Mixture (precooled, 4°C) and 50 μ l Water Phase.
- Vortex thoroughly for 5 minutes at 4°C on a vortexer with fixation aperture.
- Dispense emulsion to three thin-walled reaction tubes ("triplicates") and run reaction.

B

BINDING STEP - DNA PURIFICATION

4

4. Emulsion Breaking and DNA Binding Step

- Apply 40 μ l Activation Buffer DX to spin-column membrane, keep column at room temperature until usage, do not spin.
- Preheat Elution Buffer to 80°C (see section D, Elution).
- Pool triplicates of each sample into a 2 ml plastic reaction tube.
- Add 1 ml 2-butanol (or butanol) and break emulsion by vortexing.
- Add 400 μ l of orange-colored Orange-DX buffer (max. 250 μ l water phase / 25 μ g DNA).
- Mix buffer completely with sample.
- Centrifuge for 2 min at 11 000 x g (approx. 12 000 rpm).
- Remove organic phase (leave a small remain, do not remove interphase).
- Transfer water phase, interphase and remains of organic phase into a spin-column / receiver tube assembly.
- Centrifuge for 1 min at 11 000 x g (approx. 12 000 rpm).



C

WASHING STEPS - DNA PURIFICATION

5

5. First Washing Step

- Discard flow-through and place back spin-column.
- Add 500 μ l of Wash-DX1 buffer to spin-column.
- Centrifuge for 1 min at 11 000 x g (approx. 12 000 rpm).

6

6. Second Washing Step

- Discard flow-through and place back spin-column.
- Add 650 μ l of Wash-DX2 buffer to spin-column.
- Centrifuge for 1 min at 11 000 x g (approx. 12 000 rpm).

7

7. Removal of Wash Buffer Traces

- Discard flow-through and place back spin-column.
- Centrifuge for 2 min at 11 000 x g (approx. 12 000 rpm) to remove any remaining traces of Wash-DX2 buffer.



D

ELUTION STEP - DNA PURIFICATION

8

8. DNA Elution

- Place spin-column in new collection tube (1.5 - 2 ml).
- Add 50-150 μ l Elution-DX buffer (optional: heated to 80°C).
- Incubate for 2 min at room temperature.
- Centrifuge for 1 min at 11 000 x g (approx. 12 000 rpm).
- Discard spin-column, cap the collection tube.
- DNA is ready for analysis / manipulation or scale-up emulsion reaction.
- Store DNA at -20°C (recommended) or at 4°C (short term only).



THIS PROTOCOL IS AVAILABLE ONLINE

WWW.ROBOKLON.COM