FAST SG qPCR Master Mix (2x) + ROX REAL TIME PCR KIT

Kit Components:

Fast SG qPCR Master Mix (2x) plus ROX Solution

Component	Cat. No. E0412-01	Cat. No. E0412-02	Cat. No. E0412-03
	100 reactions, 25 µl each, 2.5 ml (1x) final volume	200 reactions, 25 µl each, 5 ml (1x) final volume	1.000 reactions, 25 µl each, 25 ml (1x) final volume
Fast SG qPCR Master Mix (2x)	1 x 1.25 ml	2 x 1.25 ml	10 x 1.25 ml
ROX Solution, 25 µM	20 µl	35 µl	130 µl
UNG (uracil-N- glycosylase) 1 U/µl	30 µl	55 µl	270 µl
Water, nuclease free	1 x 1.25 ml	2 x 1.25 ml	10 x 1.25 ml

Storage

Store at 4°C in the dark. The performance of the master mix is guaranteed up to 9 months. Freezing and thawing of the master mix should be avoided.



Figure 1: FAST SG qPCR Master Mix (blue colored curves) and "standard" SG qPCR Master Mix (green colored curves) differ by -1 Ct.
PCR conditions: 45°C 120s (UNG step) - 95°C 180 s - 45 x (96°C 5 s, 64°C 10 s, 72°C 20 s).
Ct values: Fast SG qPCR Master Mix 20.41 and 20.65, respectively, "standard" SG qPCR Master Mix: 21.73 and 21.89, respectively.

Description:

- Optimized fast qPCR master mix for 2-step and 3-step qPCR protocols.
- Fast SG qPCR Master Mix (2x) is a universal solution for quantitative real-time PCR and two-step real-time RT-PCR and is compatible with most real-time PCR cyclers available.
- The master mix contains Perpetual *Taq* DNA Polymerase, optimized reaction buffer, dNTPs (dTTP is partially replaced with dUTP) and SYBR Green I dye.
- Perpetual *Taq* DNA Polymerase contains a recombinant *Taq* DNA Polymerase bound to anti-*Taq* monoclonal antibodies that block polymerase activity at moderate temperatures.
- The polymerase activity is restored during the initial denaturation step, when amplification reactions are heated at 95°C for at least two minutes.
- Use of the "hot start" enzyme prevents extension of misprimed products and primer-dimers during reaction setup leading to higher specificity and sensitivity of PCR reactions.
- The polymerase - antibody enzyme preparation enables convenient reaction setup at room temperature.
- SYBR Green I enables DNA analysis without synthesizing expensive sequence-specific labeled probes.
- SYBR Green I is a fluorescent dye which binds all double-stranded DNA molecules and emits a fluorescent signal on binding. The excitation and emission maxima of SYBR Green I are at 494 nm and 521 nm, respectively, which are compatible with use on any real-time cycler.
- Fast SG qPCR Master Mix (2x) contains dUTP, which partially replaces dTTP. It allows the optional use of an uracil-N-glycosylase (UNG) to prevent carryover contamination between reactions. UNG removes uracil from any dU-containing contaminating amplicons, leaving abasic sites and making DNA molecules susceptible to hydrolysis during the initial denaturation step. Any introduced dUTP labelled DNA (such as dUTP labelled PCR amplicons, which may eventually accumulate in laboratories over time due to aerosol formation) is efficiently removed, while template DNA without U-residues remains unaffected by UNG treatment.
- There are two variants of the kit: without ROX and with ROX Solution provided separately. The use of ROX passive reference dye is necessary for all real-time PCR cyclers from Stratagene. ROX compensates for variations of fluorescent signal between wells due to slight differences in reaction volume and fluorescence fluctuations. ROX is not involved in PCR reaction, has a different emission spectrum than SYBR Green I and does not interfere with real-time PCR on any instrument. Refer to the table below to determine the recommended amount of ROX (25 µM) required for a specific PCR cycler.

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